


SecurAcath Subcutaneous Catheter Securement System

- Can be added onto any catheter (3F-12F)
- Small, flexible feet placed just beneath skin at insertion site
- No adhesives needed for securement
- No sutures or additional punctures needed
- Remains in place for life of catheter
- Secures catheter right where it enters the skin
- Eliminates catheter movement and pistoning
- Significantly reduces risk of catheter-related infection
- Dramatically reduces catheter dislodgement and migration
- Reduced complications mean lower total cost of patient care

Learn more at securacath.com



securAcath[®]